

# Shoprider

## 778HD

ASK FOR THE BRAND **SHOPRIDER**



**SHOPRIDER®**  
[www.shoprider.com.au](http://www.shoprider.com.au)  
*mobility*

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 778HD

**SHOPRIDER®**  
www.shoprider.com.au *mobility*

The SHOPRIDER™ 778HD is another fantastic scooter from the extensive range of scooters and powerchairs by SHOPRIDER. The 778HD is the workhorse of the three wheel range of scooters offered it has all the features that expect from Shoprider, fully rotating Airflow comfort seat,, suspension, adjustable steering tiller and the added bonus of flat free tyres as standard. This scooter also a powerful 2hp motor. The Reliability of this impressive scooter by SHOPRIDER is further enhanced by the rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

Length	1230mm
Width	640mm
Height	1250mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	182Kgs
Charger	4Amp
Motor Size	2.0Hp
Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 30km
Wheel Size	260/330mm
Turning Radius	1290mm
Weight (w/o batteries)	76kg

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

### FEATURES


High Back Seat  
Adjustable Steering Tiller  
Adjustable Width Arm Rests  
Fold Arm Rests  
Sliding Rotating Seats  
Front Shopping Basket  
Variable Seat height  
Horn  
Battery Level Warning Light  
Flat Free Tyres  
Front & Rear Suspension  
Front Lights  
Anti Tips Wheels  
**Comprehensive Warranty**

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

**WARNING:** Do not exceed 10kph on public roads or paths

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



**SHOPRIDER®**  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au