Shoptider 888SE

ASK FOR THE BRAND **SHOPRIDER** Available Colours





The SHOPRIDERTM **888SE** is packed with all the features of the larger models with the advantage of its compact size. The 888SE equipped with large batteries and has the following features, 260mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold up arm rests. The reliability of this impressive scooter by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

888SE **SPECIFICATIONS**

Length	1280mm
Width	600mm
Height	1070mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	4Amp
Motor Size	0.75Hp

Battery Capacity (Minimum Amp)	2 x 50Amp
Maximum Range #	up to 40km
Wheel Size	260mm
Turning Radius	1500mm
Weight (w/o batteries)	76kg

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia and is manufactured

to comply with all relevant international ISO standards.

IMPORTANT: Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

- # Driving range varies with the factors shown below
- (I) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

FEATURES

High Back Seat Adjustable Steering Tiller Adjustable Width Arm Rests Fold Arm Rests **Sliding Rotating Seats** Front Shopping Basket Variable Seat height Horn Battery Level Warning Light Front & Rear Bumper Pads Front & Rear Suspension Front Lights Anti Tips Wheels **Comprehensive Warranty**



*(Single user)

Your Authorised dealer



For your nearest Distributor call 1800 672 633 Or visit our Web Site at www.shoprider.com.au 6° - 113kg

0° - 136kg