

# Shoprider VENTURER

ASK FOR THE BRAND SHOPRIDER



**SHOPRIDER**<sup>®</sup>  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## SHOPRIDER<sup>®</sup> VENTURER

www.shoprider.com.au *mobility*

The SHOPRIDER™ VENTURER sets new standards for design and function for the SHOPRIDER™ range. This scooter offers a modern streamline design with many updated features including an easy to use tiller dash, full LED lighting with super bright headlight, gas adjustable tiller and 265 mm pneumatic tyres. The **Venturer** is also fitted with comfort ride suspension and a high back seat. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by a powerful 1.0hp motor, rear wheel drive and electromagnetic brake.

### SPECIFICATIONS

	VENTURER
Length	1250mm
Width	600mm
Height	1160mm
*Safe Climbing Angle	9 Degrees
*Max Climbing Angle	9 Degrees
Maximum Load Single User*	136Kgs
Charger	5Amp
Motor Size	1.0Hp
Battery Capacity (Minimum Amp)	50Amp
Maximum Range #	up to 40km
Wheel Size	265mm
Turning Radius	1440mm
Weight (with battery)	95kg

### FEATURES

- High Back Seat
- Gas Adjustable Steering Tiller
- Easy to Use Dash Panel
- LED Lighting
- Super Bright Headlight
- Variable Seat Height
- 265mm Pneumatic Tyres
- Sliding & Rotating Seat
- Adjustable Width Arm Rests
- Front & Rear Suspension
- Front Shopping Basket
- Hazard Lights / indicators
- Horn
- Anti Tips Wheels
- Comprehensive Warranty**

Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)


**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style

### VENTURER



**SHOPRIDER<sup>®</sup>**  
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au