

# Shoprider

## 889 Series

ASK FOR THE BRAND SHOPRIDER



**889SL**

**889ERS  
889ERSF**

**SHOPRIDER®**  
*mobility*  
[www.shoprider.com.au](http://www.shoprider.com.au)

RELIABLE MOBILITY GUARANTEED

# Shoprider

## 889 Series

**SHOPRIDER**  
www.shoprider.com.au mobility

The SHOPRIDER™ 889 Series provides the comfort of a large scooter with two levels of design features. Every model in the 889 Series has the same basic features, 330mm tyres for safety & comfort, fully rotating lockable seat, adjustable steering tiller and fold down arm rests. The reliability of these impressive scooters by SHOPRIDER™ is further enhanced by the rear wheel drive and electromagnetic brake.

SPECIFICATIONS	889ERS	ERSF	889SL	FEATURES	ERS	ERSF	889SL
Length	1300mm		1300mm	Large Captain's Seat	Y	Y	N
Width	580mm		640mm	Double Pillow Comfort Seat	N	N	Y
Height	1000mm		1200mm	Adjustable Steering Tiller	Y	Y	Y
*Safe Climbing Angle	12 Degrees		12 Degrees	Adjustable Width Arm Rests	Y	Y	Y
*Max Climbing Angle	15 Degrees		15 Degrees	Fold Arm Rests	Y	Y	Y
Maximum Load Single User*	136/159Kgs		136Kgs	Sliding Rotating Seats	Y	Y	Y
Charger	4Amp		4Amp	Front Shopping Basket	Y	Y	Y
Motor Size	0.75/1.5Hp		1.20Hp	Modesty Panel	N	N	Y
				Variable Seat height	Y	Y	Y
Battery Capacity (Minimum Amp)	36/50Amp		40Amp	Horn	Y	Y	Y
Maximum Range #	up to 40km		up to 40km	Battery Level Warning Light	Y	Y	Y
Wheel Size	330mm		330mm	Front & Rear Bumper Pads	Y	Y	Y
Turning Radius	1370mm		1440mm	Front & Rear Suspension	Y	Y	Y
Weight (with battery)	99kg		99kg	Front & Rear Lights	N	Y	Y
				Hazard Lights	N	Y	Y
				Front Lights	Y	Y	Y
				High / Low Speed switch	N	Y	Y
				Anti Tips Wheels	Y	Y	Y
				<b>Comprehensive Warranty</b>	<b>Y</b>	<b>Y</b>	<b>Y</b>


Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**WARNING:** Do not exceed 10kph on public roads or paths.

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

# Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



\*(Single user)

889 ERS/ SL	889 ERSF
12° (20%) 70kg	
9° (15%) 90kg	9° (15%) 113kg
6° (10%) 113kg	6° (10%) 136kg
0° (0%) 136kg	0° (0%) 159kg

**SHOPRIDER**  
www.shoprider.com.au mobility

Your Authorised dealer

For your nearest Distributor call 1800 672 633  
Or visit our Web Site at www.shoprider.com.au