



NEW PRODUCT RELEASE !



STANDARD
Bed Rail



SLOPE
Bed Rail



PEAK
Bed Rail

COMING SOON !

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

BED RAIL

RG611 - STANDARD BED RAIL

RG612 - SLOPE BED RAIL

RG614 - PEAK BED RAIL (COMING SOON)

REGISTERED DESIGNS

FEATURES

Assistive grab rail to aid user to get in and out of bed comfortably, securely and helping them to maintain their independence

Left or right hand side bed positioning
Top rail on SLOPE + PEAK bed rail can be flipped to suit

Easy to assemble - no tools required

Easy to locate and remove for bed making

Strong lightweight construction
Epoxy powder coated mild steel frame

Can pull apart for storage and transport

Optional tie down straps available (RG61BS)



SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions:

Length: **600mm**

Length: **550mm**

Length: **495mm**

Product Weight:

Approx. 2kg

Please contact
REDGUM Brand
if you require
further information
on this product.

P: 08 9248 4180

F: 08 9248 9452

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090
Part of the FOR-DE Group Pty Ltd

redgumbrand.com.au