

Presenting the Cobalt Rotating Bed

COBALT
HEALTH

Advanced support

The Cobalt Rotating Bed is designed to promote independence and provide advanced levels of support, comfort and mobility for individuals who need some assistance getting in and out of bed.

The profiling chair function sees the entire bed surface rotate smoothly, providing full-body support and minimising both user and carer strain, whilst reaching a sitting position for easier bed transfers.



Cobalt Rotating Bed

Code: B26-RB

PRODUCT FEATURES

Package includes Rotating bed, head & footboards, fold down side rails, mattress and mattress cover (sheet)

Electronic height adjustment range 35cm – 70cm (not including mattress)

Functions include profiling chair, height, backrest, knee bend, Trendelenburg and Anti-Trendelenburg

Optional 20cm bed length extension, including bolster and cover



Fabric bed end covers



Replacement cover sheet



Foot end extension



Extension bolster



Increased independence

The Rotating Bed can assist users to maintain their independence by supporting them to get in and out of bed without the need for any additional family or caregiver intervention.

Where caregiver assistance is required the Rotating Bed can help to reduce strain on the carer by electronically reaching the sitting position and making transfers easier and safer.

Bed dimensions 98 x 215cm
Mattress dimensions 90 x 190 x 8cm

Homestyle wooden beech head and footboards

Optional fabric slip on head & footboard covers in blue or grey

SWL 185kg (max user weight 165kg)

Warranty Frame 5 years, Electrics 2 Years, Mattress 1 Year

HELPING
CARE
HAPPEN