

# TREK BRAVO

## 5 Year Warranty

ENGINEERED FOR  
HIGH PERFORMANCE  
AND RESPONSIVE HANDLING

POWER, AGILITY AND A SUPERB DRIVING EXPERIENCE

The Bravo is the latest dynamic newcomer to our range. It's also our most highly refined, efficient and advanced scooter ever. Equally suited to hill climbing or shopping areas the Bravo has excellent low speed control and maneuverability.

Offering an unrivaled blend of comfort, refinement and class leading practicality that very few scooters can match. What is not matched, is our stunning 5 year warranty

Complies with new proposed Australian Standards

*NOW WITH SEAT SUSPENSION*



### Specifications\*

Length	1400MM
Width	686MM
Max Load (On Flat Ground)	185KG
Max Speed programmable up to	15KM/H
Travel range up to...	45KM
Suspension front & rear	YES
Tyre Sizes	330MM
Battery	55AH
Charger	5A
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	140A

\*subject to change without notice

### FEATURES

- \* Seat Belt
- \* Free Rear Basket
- \* High/Low Speed Switch
- \* Reclining Back Rest
- \* LED Lighting & Indicators
- \* Front and Rear Suspension

