

TREK Bravo Sport

INNOVATION PREVIOUSLY ONLY DREAMED ABOUT IS NOW A REALITY

Following in the footsteps of the proven Bravo platform, Trek has taken the ALL NEW Bravo Sport to the next level with its innovative styling, ergonomically designed seating and technically advanced features. Embodying every technical advancement Trek has conceived to date, the brand new Bravo Sport is simply Scooter Royalty. We have pushed the boundaries of what can be accomplished with precision CNC engineering to achieve what is arguably the most progressively engineered scooter ever created.

UNLESS YOU DRIVE ONE YOU WILL NOT KNOW HOW GOOD IT IS



Specifications*

Length	1445mm
Width	686mm
Max Load (On Flat Ground)	228kg
Max Speed programmable up to	15km/h
Travel range up to...	45km
Suspension front & rear	YES
Tyre Sizes	330mm
Battery	75ah
Charger	8a
Controller PG S-Drive	140a

*subject to change without notice

Features

- * Seat Belt
- * Digital Control Panel
- * Phone Holder
- * Seat Suspension
- * Lockable storage box
- * High/Low Speed Switch
- * LED Lighting & Indicators
- * Front & Rear Suspension
- * Bluetooth, USB, Speakers