


# HEIGHT ADJUSTABLE BED RAIL



## FEATURES


**RG613**

- Assistive grab rail to help the user get in and out of bed as well as sit upright
- Height adjustable
- Padded hand grip
- Easy to assemble - NO tools required
- Easy to locate and remove for bed making
- Strong lightweight construction
- Can be pulled apart for storage and transport

*Please turn over for product specifications*



SCAN THE **QR** CODE TO  
CHECK OUT THE FULL  
**REDGUMBRAND** RANGE!




Please contact **REDGUM Brand**  
if you require further information on  
this product:

P: **08 9248 4180**  
E: **orders@for-de.com.au**


**ASK FOR THE BRAND**

**redgumbrandmobility.com.au**


(All specifications may vary (+/- 10mm))




Overall width  
**600mm (W)**



Overall depth  
**485mm (W)**




Handle diameter  
**32mm (Ø)**




Handle width  
**165mm (W)**

Handle height  
6 settings - 25mm increments

- Hole 1: **545mm**
- Hole 2: **520mm**
- Hole 3: **495mm**
- Hole 4: **470mm**
- Hole 5: **445mm**
- Hole 6: **420mm**



Product Weight  
**2kg**



**ASK FOR THE BRAND**

**redgumbrandmobility.com.au**