

Protect yourself
from the sun



**CUSTOM EMBROIDERY
AVAILABLE**

**ADDITIONAL SUPPORT
STRAPS INCLUDED**

BAG CARRIER SOLD SEPERATELY

SOFT SUN CANOPY

X-SSC

FEATURES

- Provides shelter from the elements
- Specifically designed for mobility scooters
- Durable fabric canopy
- Height Adjustable
- Easy to install and remove

IMPORTANT NOTE

Requires one of the following:

- Bag Carrier (X-BC)
- Walker Carrier (X-WC)
- Side Folding Walker Carrier with Sun Canopy Mount (X-SFWCSC)

X-RANGE

P: 08 9248 4180


E: orders@for-de.com.au

ASK FOR THE BRAND

scooterextras.com.au

SPECIFICATIONS


(Dimensions are nominal and may differ +/- 10mm)



Overall width (W)
625mm

Pole width (P-W)
420mm

Embroidery area height (H)
130mm




Front width (F-W)
545mm

Overall length (L)
890mm



Product weight
3.8kg



Overall height (H)
5 settings

1165mm 1090mm
1140mm 1065mm
1115mm

Measured from the scooter socket (eg. where the Bag Carrier inserts) to the underside of the back of the X-SSC.

Bag Carrier sold separately.