

Lightweight Folding Powerchair

RC-Inn-N5909



ROTHCARE



Features

- Super Lightweight carbon fibre frame
- Folds in one easy step for transportation
- Comfortable padded seat & backrest
- Lightweight removable Lithium battery
- Compact easy to use joystick
- Complete chair weighs less than 13kg with the battery
- USB fitting in Joystick & in removable battery pack
- Battery pack is fitted with a charging port for offboard charging

Tech Specs

| | |
|----------------------|---------------------------|
| Motor | DC 24V, 150W brushless DC |
| Controller | Shanghai Micon |
| Brake | Electromagnetic |
| Front wheel | 6" |
| Rear wheel | 10.6" |
| Weight capacity | 100kg |
| Weight w/out battery | 11.85kg |
| Weight with battery | 13kg |

| | |
|-------------------------|------------------------|
| Overall width | 53cm |
| Overall length | 90cm |
| Overall height | 86cm |
| Seat width | 43cm |
| Seat depth | 40cm |
| Seat height from ground | 46cm |
| Backrest height | 40cm |
| Speed | 4.5kph |
| Turning radius | 90cm |
| Gradient | 3 degrees |
| Range | 10km |
| Battery | 24V / 10Ah Lithium-ion |