


# DRESSING STICK

RG7IDS

## METAL ROUND HOOK

that can be hung when not in use



Assist with  
**DRESSING AND  
REMOVING**  
clothing

Solid wood  
handle

Dual prong  
rubber hook

Dual prong  
rubber hook

## FEATURES




- Assist with dressing and removing clothing
- Solid wood handle
- Dual prong rubber hook
- Metal round hook
- Can be hung when not being used

## SPECIFICATIONS

- Length: **690mm**
- Handle: **Ø15mm**
- Hooks: **Ø6mm / Ø3mm**
- Product Weight: **70g**

Please contact **REDGUM Brand** if you require further information on this product.

**P: 08 9248 4180**  
**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**



Hook can be  
used to zip  
up jacket

**ASK FOR THE BRAND**

**[redgumbrand.com.au](http://redgumbrand.com.au)**