



Eco 745

 **merits**
AUSTRALIA

www.meritsaustralia.com


Eco 745

Features

- Stylish eye catching design with attractive modern colour range
- Loop handlebars offer comfortable finger or thumb throttle operation
- Adjustable high back contoured seat with simple operation of slide and swivel mechanisms to allow maximum user comfort
- Infinite adjustment of steering tiller allowing everyone to easily find a comfortable position
- High mount charger socket means easy access when charging
- Removable front basket with carry handle
- Adjustable power saving LED headlight, indicators front and rear, rear tail lights and reverse warning
- Built in speed reduction switch which is activated when turning offering exceptional safety.



Dimensions cm/in



Overall Dimension A x B x C	120 x 61 x 99cm/47 x 24 x 39"
Seat Dimension : D x E x F	46 x 40 x 46cm/18 x 16 x 18"
Seat Height from Ground : G	53-58cm/21-23"
Seat Height from Deck : H	40-47cm/16-18.5"
Wheel Base : I	85cm/33.5"

Weight kg/lb

base	battery/ea	seat	total with battery	capacity
51/112	10/22	15/33	86/190	136/300

Technical Specs

Max Speed	10kph/6.25mph
Range	20-30km/13-20mi
Turn Radius	157cm/62"
Ground clearance	9.5cm/3.74"
Motor	DC 24V, 210W (Rating)/715W (Max)
Controller	S-Drive 70A
Battery	12V/34A/H x 2pcs
Charger	5A off-board
Gradient	6°
Front Wheel	10" Pneumatic tire
Rear Wheel	10" Pneumatic tire
Brake	intelligent, regenerative, electromagnetic brakes

merits
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.