

# TREK Evolution 4

COMBINING ELEGANT SIMPLICITY WITH AGILITY  
AND UNMATCHED RELIABILITY

UNRIVALLED COMFORT REFINEMENT AND PERFORMANCE

Possibly the best overall luxury scooter value in its class. The Evolution 4 is ideally suited to users who require a feature rich luxury product at an attractive price.

Designed by Korea's leading scooter manufacturer, features include a powerful 4 pole motor, front and rear suspension, LED lighting and fully adjustable high back seat.

**Phone Holder - Bluetooth - Speakers - USB**

*-Charge your phone*

*-Play music via Bluetooth or USB stick*

*-Stream radio stations from your phone*



## Specifications\*

Length	1345MM
Width	695MM
Max Load	182KG
Max Speed programmable up to	12KM/H
Travel range up to...	40KM
Suspension front & rear	YES
Tyre Sizes	300MM
Battery	45AH
Charger	6A
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	120A

## Features

- \* Seat Belt
- \* Free Rear Basket
- \* High/Low Speed Switch
- \* Reclining Back Rest
- \* LED Lighting & Indicators
- \* Front and Rear Suspension
- \* Seat Suspension

\*subject to change without notice

trekmobility.com.au