

TREK Evolution

COMFORT & RELIABILITY AT AN AFFORDABLE PRICE

Simple designs and high quality components translate to a long life cycle and very low running costs. Our high quality Taiwanese drive system has a proven track history and will provide you with years of maintenance free low cost and worry free travel. A high level of comfort is guaranteed for those long trips regardless of surface conditions with our newly designed seat suspension system and high back reclining seat.

Phone Holder - Bluetooth - Speakers - USB

- Charge your phone
- Play music via Bluetooth or USB stick
- Stream radio stations from your phone



Specifications*

Length	1270mm
Width	620mm
Max Load (On Flat Ground)	160kg
Max Speed programmable up to	10km/h
Travel range up to...	35km
Tyre Sizes	260mm
Battery	36ah
Charger	6a
Motor (Taiwanese 2 Pole)	
Controller PG S-Drive	120a

Features

- * Seat Belt
- * Free Rear Basket
- * High/Low Speed Switch
- * Reclining Back Rest
- * LED Lighting & Indicators
- * Seat Suspension

*subject to change without notice

trekmobility.com.au