

SCOOTER XTRAS

IMPROVED
DESIGN !



Mobility Scooter
not included



Bag + Walking Stick Carrier (SE-BC)
required to fit Hard Canopy
(Sold separately)

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

HARD CANOPY

SE-HC

FEATURES

- Provides excellent wind + sun protection
- Clear polycarbonate windscreens
- Lightweight metal framework
- Easily removable for servicing

Hard Canopy can only be fitted to the following
LISTED SHOPRIDER mobility scooters:
(2012 models and later only)

ROCKY 8 (889AE) **SE-HCMKR8 Mount kit required**

ROCKY 6 (889XLSBN + B) Mount kit included

ROCKY 4 (889XLB) Mount kit included

ALLROUNDER (9LSP) **SE-HCMK2 Mount kit required**

889SL (889SL) **SE-HCMK2 Mount kit also required**

Bag + Walking Sticker Carrier (SE-BC)
or

Wheelie Walker Carrier (SE-WWB)

is required to fit Hard Canopy to any scooter
(Sold separately)

Optional Accessories (Sold separately)

Hard Canopy Side Covers (SE-HSSC)

- Clear PVC side and rear enclosure curtains

SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions:

1300mm (L) x 720mm (W) x 1400mm (H)

Product Weight: **8kg**

NOT SUITABLE ON ANY SCOOTER PRE 2012

Warning:

Using the Hard Canopy in strong winds and adverse weather can seriously compromise your safety and may damage the mobility scooter beyond repair. The hard canopy is NOT to be used on a scooter on a GOLF COURSE situation as damage to the hard canopy and scooter will occur.

Please contact
 Scooter Xtras
if you require
further information
on this product



Hard Canopy Side Covers (SE-HSSC)
(Sold separately)

P: 08 9248 4180

E: orders@for-de.com.au

Scooter Xtras - 1 Business way, Malaga, Western Australia 6090
Part of the **FOR-DE Group** Pty Ltd

scooterxtras.com.au