

# 4PT™

4 POLE TECHNOLOGY



S940  
Interceptor  
Series 2

merits

# S940 Interceptor Series 2



## Features

- Brand new sporty streamline design.
- Fully adjustable suspension system that delivers a smooth and comfortable ride.
- Modern LED lighting system delivers a brilliant light with less power consumption.
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer.
- Extra storage space beneath tiller.
- Built-in half speed switch which operates automatically when turning to provide extra safety.
- With Merits 4PT (4 pole technology) and large controller delivering much more power for climbing steep driveways and hills.
- With new easy to access 12 volt socket.



Modern LED lighting system delivers brilliant light with less power consumption.




More space for storage


Digital LCD dashboard display

Handy 12 volt socket

## Dimensions cm




seat width  
51cm




seat depth  
46cm



seat height  
from ground  
64-71cm



seat height  
from deck  
47-52cm



back height  
63cm



wheel base  
113cm



length  
154cm



width  
72cm



height  
120cm

## Weight kg

base	battery/ea	seat	total with battery	capacity
84	25	18	165	204

## Technical Specs.

speed	10-14kph
range up to	50km
turn radius	180cm
ground clearance	13cm
motor	DC24V,2200W, 4 pole
controller	PG S-Drive 140A
battery	12V / 60-75AH x 2
charger	8A off-board
gradient	10°
front wheel	14" (3.50-8) pneumatic tyre
rear wheel	16" (4.00-10) pneumatic tyre
brake	intelligent, regenerative, electromagnetic brakes.

Scooter World Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.



[www.scooterworldaustralia.com.au](http://www.scooterworldaustralia.com.au)