

PRIDE[®]
MOBILITY



<20 km
Range per Charge



136 kg
Weight Capacity



56 cm
Turning Radius



Jazzy[®]
EVO 613

The Most Advanced Jazzy[®] Ever

The advanced comfort and versatility of the Jazzy[®] EVO 613 provides the freedom to navigate small spaces and tight corners. Its narrow base width of only 56cm makes it an ideal powerchair for moving through your home. The high-back swivel seat provides convenience and easy access, with the additional option of a swing-away controller mount



Visit pridemobility.com.au for more information

Jazzy[®] EVO 613



Narrow Base Width of 55.8 cm



Swing Away Controller Mount (optional extra)



Mid-Wheel 6[®] Technology

| | |
|-------------------------------------|---|
| RANGE PER CHARGE (MAX) ¹ | Up to 20 km |
| MAXIMUM WEIGHT CAPACITY | 136 kg |
| TURNING RADIUS | 56 cm |
| MAXIMUM SPEED ¹ | 7.8 km/h |
| LENGTH & WIDTH ² | L: 104 cm W: 55.8 cm |
| WEIGHT OF HEAVIEST PIECE | 42.6 kg |
| GROUND CLEARANCE | 7.6 cm at motors |
| SEAT SIZE WIDTH & DEPTH | W: 45.7 cm x D: 45.7 cm (18" x 18" - 20") |
| STANDARD 18" X 18" SEAT WEIGHT | 22.7 kg |
| HIGH-BACK SEAT-TO-FLOOR RANGE | 50.2 cm - 54 cm |

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Not suitable for air travel

| | |
|-----------------------------------|--|
| STANDARD ELECTRONICS | 40A, PG GC3 |
| DRIVE TRAIN | Two motor, in-line, mid-wheel drive |
| BRAKING SYSTEM | Regenerative and electro-mechanical |
| SUSPENSION TYPE | Active-Trac [®] (ATX) |
| FRONT CASTORS | 5" solid |
| DRIVE WHEELS | 13" solid |
| REAR CASTORS | 5" solid |
| BATTERY REQUIREMENTS ³ | (1) (2) 31Ah, 12V |
| BATTERY CHARGER | Off-board, 3.5A |
| BATTERY WEIGHT | 8.4 kg |
| COLOUR OPTIONS | <div style="display: flex; align-items: center; gap: 10px;"> <div style="display: flex; align-items: center; gap: 5px;"> Black </div> <div style="display: flex; align-items: center; gap: 5px;"> Robins Egg Blue </div> <div style="display: flex; align-items: center; gap: 5px;"> Cardinal Red </div> </div> |

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.

Transport of this item on commercial aircraft is **not permitted** for the Jazzy Evo 613 Li.