

# Knee Walker

## RC-HCT-9125-F

.....



# ROTHCARE

### Features

.....

- Stylish folding frame with folding tiller for compact transportation
- Large pneumatic wheels for easy use indoors or outdoors
- Handy storage basket & comfortable knee cushion
- Handle height & knee cushion height adjustment
- Soft touch breaks with parking brake function
- Tyre pump included

### Tech Specs

.....

Length	90cm
Overall width	44cm
Knee Cushion height	47cm-60cm
Handle height	86cm-107cm
Knee Cushion thickness	7cm
Weight	11.7kg
Weight Capacity	160kg
Front Wheel (Pneumatic)	12½" x 2¼"
Rear Wheel (Pneumatic)	10" x 2"

