

INCLUDES WIRE BASKET!

Height adjustable
contour padded
kneepad



KNEE WALKER

RG40KWHD

FEATURES

- Height adjustable contour padded kneepad
- Height adjustable handlebars
- Ergonomic rubber hand grips
- Large wheels for excellent mobility
- Dual hand brakes
- Easy to assemble
- Strong lightweight construction
- Removable wire basket
- Walking stick holders
- Folds down for transport & storage

SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

Overall:

500mm (W) x 1050mm (L) x 900mm (H)

Kneepad Size:

190mm (W) x 380mm (L)

Kneepad Height: (from floor to top of seat)

530mm - 640mm

Handlebar Width: **490mm**

Handlebar Height: (from floor to top of handlebar)

900mm - 1050mm

Wheel Diameter: **200mm / 8"**

Product Weight: **15.5kg**

Max. User Weight: **180kg**



Folds down
for transport
& storage

Please contact **REDGUM Brand**
if you require further information
on this product.

P: 08 9248 4180

E: orders@for-de.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

redgumbrand.com.au