



**NEW
PRODUCT
RELEASE!**



LIFT
position



RECLINE
position



TOFFEE
(LC0301)



GREY
(LC0302)



BLACK
(LC0303)

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

FLORENCE

WALL SAVER LIFT CHAIR

LC0301 - TOFFEE | LC0302 - GREY | LC0303 - BLACK

FEATURES

Compliant to AS/NZ 4688

- Lift / Rest / Recline armchair.
- Allows the user to be easily raised to a standing position and offers the comfort of a backrest recline.
- Chair can be positioned close to wall (150mm).
- Single motor lift.
- Luxurious leather coverings.
- Available in modern colours: Toffee / Grey / Black.
- Quality German actuator with battery backup.
- Warranty:
Frame: 3 years / Electrics: 2 years / Fabric: 2 years

SPECIFICATIONS

All dimensions shown are nominal

Overall: **740mm (W) x 970mm (D) x 1020mm (H)**
Seat: **480mm (W) x 480mm (D) x 480mm (H)**
Armrest: **150mm (H) x 460mm (L)**
Backrest: **670mm (W) x 650mm (H): 93 - 125°**
Legrest: **540mm (W) x 480mm (L): 90 - 180°**
Distance from wall: **150mm**
Product Weight: **52kg**
Maximum User Weight: **160kg**

Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180
E: orders@for-de.com.au

1 Business way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au