

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

TARANTO FABRIC 4 MOTOR LIFT CHAIR

LC0401 - DARK GREY

LC0402 - TAUPE

FEATURES

- Lift / Rest / Recline armchair.
- Independent headrest and lumbar adjustment.
- Allows the user to be easily raised to a standing position and offers the comfort of a full recline.
- Four independent motors allow you to maximise vour comfort.
- Added comfort of a pocket spring seat base.
- Stain resistant fabric coverings: Dark Grey / Taupe.
- Quality German actuators with battery backup.
- Warranty:

Frame: 5 years / Electrics: 2 years / Fabric: 2 years

SPECIFICATIONS All dimensions shown are nominal

Overall: 905mm (W) x 940mm (D) x 1080mm (H) Seat: 470mm (W) x 500mm (D) x 500mm (H) Armrest: 210mm (W) x 500mm (L) x 170mm (H) Backrest: 560mm int. (W) x 670mm (H): 93 - 180° Legrest: 620mm (W) x 430mm (L): 90 - 180° Headrest: 560mm (W) x 260mm (H): 0 - 45°

Distance from wall: 650mm Product Weight: 65kg

Maximum User Weight: 160kg

Please contact **REDGUM Brand** if you require further information

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au