

## Shower / Toilet Transport Commode

RG8923




### Features

- Powder coated aluminium frame
- Pearl white finish
- Can be used to fit over a toilet
- Secure fitting bucket with lid
- Perfect for shower usage
- Padded seat included
- Moveable armrest for easy side transfer
- Padded armrest
- Easy move 120mm rubber castors
- Lockable rear castors fitted
- Leg rests fitted - Flip back for storage
- Easy remove seat for cleaning

### Dimensions

- Seat Width 18" (450mm)
- Seat Depth 440mm
- Seat Height 530mm
- Seat Back Height 420mm
- Max User Capacity - 130kg
- Unit Weight - 10.5kg



Phone: 08 9248 4180  
Fax: 08 9248 9452  
E-mail: [orders@shoprider.com.au](mailto:orders@shoprider.com.au)