

# Pathrider® DX Sport

HIGH-PERFORMANCE OPERATION  
AT AN AFFORDABLE PRICE

## PATHRIDER® DX SPORT

The all-new compact **Pathrider® Deluxe Sport** Scooter features a sleek, sporty design and delivers high-performance operation. Just some of the standard features include a handbrake, front and rear suspension, feather-touch disassembly, a wraparound delta tiller, high-visibility automatic brake lights and an LED kerb light. These features combine to make the Pathrider DX Sport a superb scooter at an exceptional price.



**Pride**<sup>®</sup>  
Mobility Products  
Australia Pty. Ltd.

— LIVE YOUR BEST<sup>®</sup> —

(Tel) 03 8770 9600 • (Fax) 03 9703 2720 • [www.pridemobility.com.au](http://www.pridemobility.com.au)

# Pathrider® DX Sport

Weight capacity:	158.75 kg (350 lbs.)
Maximum speed <sup>1</sup> :	Up to 12.87 km/h (8 mph)
Ground clearance <sup>2</sup> :	8.25 cm (3.25")
Turning radius:	133.35 cm (52.5")
Overall length <sup>2</sup> :	119.38 cm (47")
Overall width <sup>2</sup> :	66.37 cm (26.13")
Front tyres:	9 cm x 25.4 cm (3.5" x 10")
Rear tyres:	10 cm x 25.4 cm (4" x 10")
Suspension:	Front and rear
Range per charge <sup>1,5</sup> :	Up to 16 miles (25.96 km)
Maximum obstacle climbing ability:	5 cm (2")
Total weight without batteries:	68 kg (150 lbs.)
Heaviest piece when disassembled:	Front section: 29.48 kg (65 lbs.)
Standard seating:	Type: CRS Junior seat Weight: 19.05 kg (42 lbs.) Material: Black Vinyl Dimensions: width – 44.45 cm (17.5") overall width depth – 43.18 cm (17") effective depth
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: 35-40 AH Weight <sup>4</sup> : <b>35-40 AH</b> : 14.74 kg (32.5 lbs.) each
Battery charger:	5-amp, off-board (supplied with 35-40 AH batteries)

Colours:



<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition.

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>AGM or gel cell type required.

<sup>4</sup>Battery weight may vary based on manufacturer.

<sup>5</sup>Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

## Features

- Front and rear suspension
- High-visibility automatic rear brake lighting system
- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Long-lasting, high-intensity LED headlight, taillight and indicators
- Exclusive stylish, lightweight, non-marking, low-profile wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- New LED battery meter for increased accuracy and visibility in any setting
- Includes pinchless hinges & eight-position seat rotation
- Easy-grip tiller adjustment knob
- Easy access, tiller-mounted charger port for standard dual voltage off-board charger
- Rearview mirror standard
- Handbrake

## Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder