



QUICKIE
Live without limits.

Q-Care
Powerchairs

Q300 M CAPTAIN SEAT

Compact in size, BIG on performance.

Combining over 40 years of innovative power wheelchair design



GO WHERE OTHER POWER CHAIRS CAN'T WITH THE Q300 M MINI CAPTAIN SEAT

Ultra-slim at just 520 mm, the all-new Q300 M Mini is the narrowest TRUE mid-wheel drive powerchair - and with a turning circle of just 1120 mm. This chair will take you places that other chairs can't fit.

CRASH TESTED

16" 18" & 20" WIDE

30° TILT OR 12" LIFT

COMES WITH ULTRA-SLIM BASE

Ultra-narrow 520 mm width
Ultra-compact 1120 mm turning circle

GIVES YOU SPACE WHEN YOU NEED IT

Low seat height from 470 mm
Easy transfers with flip-back armrests
Get tighter turns with the Centre Mount legrest

PERFECT FOR INDOOR NAVIGATION

12.5" Standard Drive wheels for the tightest manoeuvrability and indoor practicality

GREAT TRAVELLING COMPANION

Crash-test approved with Four-Point tie down restraint systems
Fits under the steering wheel of your wheelchair adapted vehicles

TRUE OUTDOOR PERFORMER

Narrow enough to be nimble in a crowded city
Anti-Pitch Technology keeps you safe and stable
All-wheel independent suspension for a smooth ride
Powered by punchy 50 Ah batteries, it takes you everywhere you want to go

CAPTAIN SEAT: COMFORTABLE CAR STYLE SEAT

Comfortable car style seat that provides a relaxed sitting position
Features manual recline

Seat cushion can be replaced with pressure cushion

POWER SEAT OPTIONS

Relax with 30° of power tilt
Access your environment with 12" of lift
A reliable VR2 control and a simple to control joystick which allows for intuitive operation



TECHNICAL DATA

| | | | |
|------------------|-----------------------------------|------------------------------|--|
| Max. user weight | 136 kg | Drive Type | TRUE Mid-Wheel Drive |
| Seat width | 16", 18", 20" | Max. kerb climbing | 70 mm (with 12" drive wheels) |
| Seat depth | 18" - 20" | Wheel options | 12" |
| Seat height | 18.5" - 19" | Max. range (ISO 7176-4) | 30 km |
| Backrest height | 25" | Wheelchair weight | From 123 kg (including 50 Amp) |
| Backrest angle | 90° - 115° | Crash tested | Yes - approved for occupant forward facing transport in a vehicle (ISO 7176-19). |
| Total width | 520 mm (20") | Intended use and environment | For people with limited mobility, indoor and outdoor use. |
| Total length | 1060 mm (41.5") - 1180 mm (46.5") | | |
| Tilt angle | 0° - 30° (Power Tilt) | | |
| Turning diameter | 1120 mm (44") | | |
| Max Speed | 10 kph | | |
| Battery size | 50 Ah | | |
| Electronics | VR2 | | |

SKU: Q300CAPT - SELECT OPTIONS FOR YOUR Q300 HERE

| | | |
|--|---|---|
| Seat Size <input type="checkbox"/> 16" x 16" <input type="checkbox"/> 18" x 18" <input type="checkbox"/> 20" x 18" | Leg Rests <input type="checkbox"/> Swing Away <input type="checkbox"/> Centre Mount Platform | Joystick <input type="checkbox"/> Left <input type="checkbox"/> Right |
| Power Seat Functions <input type="checkbox"/> Tilt or <input type="checkbox"/> Elevate | Shroud Colours <input type="checkbox"/> White <input type="checkbox"/> Black <input type="checkbox"/> <input type="checkbox"/> Red <input type="checkbox"/> Blue <input type="checkbox"/> <input type="checkbox"/> Purple <input type="checkbox"/> Green <input type="checkbox"/> <input type="checkbox"/> Pink <input type="checkbox"/> Silver <input type="checkbox"/> | Accent Colours <input type="checkbox"/> Green <input type="checkbox"/> Red <input type="checkbox"/> <input type="checkbox"/> White <input type="checkbox"/> Blue <input type="checkbox"/> <input type="checkbox"/> Black <input type="checkbox"/> Silver <input type="checkbox"/> <input type="checkbox"/> Orange <input type="checkbox"/> |