## Shoppider PUMA ASK FOR THE BRAND SHOPRIDER

## <image>

**PUMA 10** 

PUMA 14HD





The SHOPRIDER<sup>™</sup> PUMA has been ergonomically engineered to provide a robust combination of style and comfort. The PUMA is fitted with powerful smooth motors, mid wheel manoeuvrability, and comfortable seating. The smooth operation of the front and rear elastometers combined with mid wheel drive means easy operation on a variety of terrains and hard surfaces. It's ideal for both inside and outside use on even ground and can travel long distances between recharging.


SPECIFICATIONS	PUMA10	PUMA14HD	FEATURES
Length Width Height Battery Capacity #Maximum Range Maximum Load Single User* Charger on board	1060mm 620mm 1160mm 2 x 36Amp up to 20km 136Kgs 4Amp	1060mm 640mm 1180mm 2 x 75Amp up to 25km 205Kgs 8 Amp	Large Captain's Seat 20" Seat PUMA 14HD Sleek body design Smooth operation Variable Seat height Controller interchangeable both arms Free wheel device Adjustable arm rests Option elevating leg rests Option of swing away leg rest Stamp supports <b>Comprehensive Warranty</b>
Motor Size Maximum Speed *Safe Climbing Angle *Max Climbing Angle Turning Radius	0.8 Hp 8Kph 8 Degrees 12 Degrees 570mm	I.34 Hp 6Kph 8 Degrees 9 Degrees 600mm	
Wheel Size: - Front Caster - Middle (drive) - Rear caster	2 x 175mm 2 x 260mm 2 x 175mm	2 x 175mm 2 x 350mm 2 x 175mm	Shoprider is listed and complies with the regulations of the Therapeutic Goods Association of Australia as well as the regulatory bodies and standards in America (FDA) and Europe (CE)

**IMPORTANT:** Always be sure to read your owner's manual thoroughly before operating your machine.

WARNING: Do not exceed 10kph on public roads or paths

- # Driving range varies with the factors shown below
- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style





Your Authorised dealer

For your nearest Distributor call 1800 672 633 Or visit our Web Site at www.shoprider.com.au