

Power Pack for manual wheelchair P001F

Features

Easy attachment to manual wheelchairs to allow motorised attendant assistance

Convenient forward & reverse direction switch & speed control, with easy to operate thumb throttle

Portable & lightweight for easy transportation which enables extra mobility anywhere

Twin wheel design which allows much greater traction on grass & gravel



Technical Specs

| | |
|------------------------------|--|
| Motor | 280watt |
| Battery | 12V - 20a/h |
| Charger | 3A off board |
| Drive Wheel | 8" PU tyre x 2 pieces |
| Weight of battery | 6.3kg |
| Weight of base | 7.7kg |
| Total weight | 14kg |
| Weight Capacity | 136kg at (6 degree gradient) – 159kg at (0 degree gradient) |
| Range up to | 15km |
| Fits wheelchair width | 16inch – 22inch |

Universal mount kit
available as an option