



**FULL
HANDGRIP
TRIGGER**

**SOFT RUBBER ON JAWS
HOOK AND MAGNET AT THE TOP
IN ORDER TO MAKE IT EASIER TO
PICK THINGS UP!**

FEATURES

- Excellent aid to extend user's reach
- Reduces discomfort when bending to retrieve items
- Minimal handgrip required to operate jaws
- Comfortable squeeze lever handle action activates jaw grips
- Full hand grip trigger
- Rubber grip on jaws
- Magnet on top of jaw
- Lightweight aluminium arm with plastic parts

SPECIFICATIONS

Overall Dimensions:

RGM24:

L: 610mm x W: 40mm x D: 130mm

RGM30:

L: 760mm x W: 40mm x D: 130mm

RGM35:

L: 890mm x W: 40mm x D: 130mm

Jaw Width: **80mm**

Product weight:

RGM24: **170g** | RGM30: **185g** | RGM35: **200g**

Max. load weight: **1.5kg**

P: 08 9248 4180

E: orders@for-de.com.au

1 Business way Malaga Western Australia 6090

Part of the **For-de Group Pty Ltd**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.