



**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use.  
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.



**FEATURES**

- High quality walking stick
- Aluminium construction
- Height adjustable
- Ergonomic handle
- Convenient wrist strap
- Non slip rubber tip (37mm)
- Easy to fold for transport and storage
- No tools required

*Please turn over for product specifications*



Scan the QR code to check out our *walking stick* range!

Please contact **REDGUM Brand** if you require further information on this product:  
P: 08 9248 4180  
E: [orders@for-de.com.au](mailto:orders@for-de.com.au)




**ASK FOR THE BRAND**

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)

## Unfolded

Overall Width  
Handle - **140mm**




Overall Height  
**5 Settings (25mm)**

**850mm**  
**875mm**  
**900mm**  
**925mm**  
**950mm**

Overall Width  
Rubber tip - **37mm**

## Weight




Product Weight  
**210g**



Maximum load weight  
**113kg**

## Folded


Overall Width  
**140mm**




Overall Height  
**310mm**

## Handle

Handle Width  
**25mm**



Handle Depth  
**32mm**



Handle Height  
**60mm**

**ASK FOR THE BRAND**

[redgumbrandmobility.com.au](http://redgumbrandmobility.com.au)