

SOLD IN PAIRS!



**APPROPRIATE
FOR INDOORS
AND OUTDOORS**

Ø25mm
Internal Fit



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

SKI GLIDE

RG25 | RG25B

FEATURES

- Smoother travel over most surfaces
- Made from durable plastic
- Appropriate for indoors and outdoors
- Tool free instalation
- Smooth bottom surface
- Sold in pairs

SPECIFICATIONS


RG25

To suit: Ø28mm shaft
Footprint: 85mm (L) x 40mm (W)
Max. Load: 130kg
Colours: **Grey**

RG25B

To suit: Ø25mm shaft
Footprint: 75mm (L) x 30mm (W)
Max. Load: 130kg
Colours: **White**

Ø28mm
External fit



Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

redgumbrand.com.au