

SOLAX MOBIE

ONE STEP MANUAL FOLDING MOBILITY SCOOTER



LITHIUM BATTERY



PUNCTURE PROOF
TYRES



LIGHTWEIGHT
COMPACT DESIGN



SOLAX MOBIE

ONE STEP MANUAL FOLDING MOBILITY SCOOTER

The Mobie is the signature product of the Solax brand and what made Solax an instant success and leader in the mobility market. Although the design has been improved since it's release in 1993, the main design principles such as ease of use, safety, portability and quality are stronger than ever. An easy under-seat folding mechanism allows for easy folding and unfolding. The padded seat, adjustable tiller and flip-up removable armrests make the Mobie a truly fun and comfortable scooter.

SPECIFICATIONS

Weight With Batteries	27.3kg	Tyres	Solid - Puncture Proof
Battery Weight	1.7kg	Dimensions LxWxH (mm)	990x540x840
Heaviest Part	25.6kg	Folded Dimensions (mm)	600x540x470
Weight Capacity	125kg	Distance Between Armrests	450mm
Adjustable Tiller Angle	Yes	Drive Wheel	170mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	40mm
Battery	24V/10ah Lithium	Seat Dimensions	390mm x 340mm
Battery Range	Up To 15km	Seat Height From Ground	560mm
Maximum Speed	6km/h	Slope Ability	12°
Motor	120W	Turning Circle	1250mm
Number of Components	1 + Battery	Suspension	No

Colours Subject to Availability



Dealer Details