

SOLAX CHARGE

AUTOMATIC FOLDING MOBILITY SCOOTER



ADJUSTABLE
REAR SUSPENSION



REMOTE CONTROL
AUTOMATIC FOLDING



USB PORT FOR
CHARGING PHONE




SOLAX CHARGE

AUTOMATIC FOLDING MOBILITY SCOOTER

The Solax Charge is the first Solax that features a USB charging port that allows you to charge your phone or iPad. Other tech features include automatic folding, LED headlights, rear suspension, a digital dashboard that displays your speed, trip distance and battery life. A re-designed seat boasts a higher backrest with extra padding for additional comfort, while the armrests fold up for easier access.

SPECIFICATIONS

Weight With Batteries	26kg	Tyres	Solid - Puncture Proof
Battery Weight	1.7kg	Dimensions LxWxH (mm)	1000x550x865
Heaviest Part	24.3kg	Folded Dimensions (mm)	410x550x690
Weight Capacity	125kg	Distance Between Armrests	480mm
Adjustable Tiller Angle	Yes	Drive Wheel	170mm
Anti-Tip Mechanism	Rear Anti-Tip Wheels	Ground Clearance	38mm
Battery	24V/10ah Lithium	Seat Dimensions	430mm x 340mm
Battery Range	Up To 15km	Seat Height From Ground	530mm
Maximum Speed	6km/h	Slope Ability	12°
Motor	120W	Turning Circle	1250mm
Number of Components	1 + Battery	Suspension	Rear Adjustable



Colours Subject to Availability



DEALER DETAILS