



**FOLDS FOR
STORAGE!**



**STRONG
LIGHTWEIGHT
ALUMINIUM
FRAME**



**For more information
about this product**

FOLDING SHOWER CHAIR

RG9400

FEATURES

- Easy to assemble - No tools required
- Adjustable in height with Sure Lock Adjustment
- Strong lightweight construction - Aluminium frame
- Rubber leg tips for excellent stability in wet areas
- Reinforced plastic seat
- Flat pack - Excellent for transport + storage
- Contoured plastic seat with draining holes
- Lifetime warranty against corrosion
- Cutout plastic seat for ease of cleaning

SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

Overall: (at lowest setting - 4 settings)

530mm (W) x 500mm (D) x 880mm (H)

Between Armrests: **475mm**

Width: **410mm**

Depth: **400mm**

Height: **450mm - 525mm** (4 settings)

Backrest: **490mm (W) x 250mm (H)** - from seat

Product Weight: **4.5kg**

Max. User Weight: **100kg**

P: 08 9248 4180
E: orders@for-de.com.au

ASK FOR THE BRAND

redgumbrandmobility.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.