

ROUND BATH STOOL

RG9012 (STANDARD) | RG9012R (ROTATING)

FEATURES

- Strong lightweight construction
- Designed to make showering and cleaning easier for user
- Padded seat
- Height adjustable
- Strong aluminium legs
- Easy to assemble - no tools required

SPECIFICATIONS

(All dimensions shown are nominal and at lowest setting)

RG9012 - Standard Round Bath Stool:

- Overall: **325mm (W) x 325mm (D)**
490 - 660mm (H) (7 increments)
- Seat Diameter: **Ø325mm**
- Product Weight: **3.1kg**
- Max. Supportive Load: **150kg**

RG9012R - Rotating Round Bath Stool:

- Overall: **325mm (W) x 325mm (D)**
500 - 670mm (H) (7 increments)
- Seat Diameter: **Ø325mm**
- Product Weight: **3.5kg**
- Max. Supportive Load: **150kg**

Please contact **REDGUM Brand** if you require further information on this product

P: **08 9248 4180**

E: **orders@for-de.com.au**



**ROUND SHOWER
STOOL**
(RG9012)



**ROTATING ROUND
SHOWER STOOL**
(RG9012R)



**PADDED
SWIVEL SEAT**
for rotating round
shower stool

Locks every 90°



ASK FOR THE BRAND

redgumbrandmobility.com.au