

TREK Sun Rider



REDEFINING WHAT A MOBILITY SCOOTER
SHOULD BE. POWER, AGILITY, COMFORT
AND A DYNAMIC DRIVING EXPERIENCE

**5 Year
Warranty**

NOW WITH SEAT SUSPENSION



Shaping The Future Of Mobility

With its carefully crafted protective canopy designed to optimise usage in all weather conditions the Sun Rider stands out from the crowd. Its assertive design is complimented with a refined blend of performance and luxury making it our most advanced, efficient and refined scooter ever.

Powerful, agile and distinctive, the TREK Sun Rider is engineered for safety, high performance and responsive handling.

Complies with new proposed Australian standards

TREK SUN RIDER


Australia's Best Scooter Warranty

5 Years on Motor & Transaxle



*Complies with proposed new Australian standards
(Legal in all states including Queensland)*

Easy open access for Trouble Free Maintenance



Features

- *Reclining Back Rest
- *Glass Windscreen
- *High & Low Speed
- *Front & Rear Suspension
- *Large PU Front Bumper
- *LED Lighting
- *Indicators
- *Twin Mirrors
- *Rear Basket
- *Seat Belt



Specifications*

Length	1430MM
Width	686MM
Height	1600MM
Weight	145KG
Max Load (On Flat Ground)	185KG
Max Speed programmable up to	14KM/H
Travel range up to...	40KM
Suspension front & rear	YES
Tyre Sizes	330MM
Battery	55AH
Charger	5A
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	140A

*subject to change without notice

trekmobility.com.au