

TREK SunRider Rio

A PREMIUM SCOOTER AT AN AFFORDABLE PRICE

Inspired by the success of our flagship SunRider, our design team at TREK are proud to introduce the affordable all weather mid-range SunRider Rio to our marque range of products.

Designed to TREK's specifications, the Rio's superior torque and handling outperforms all competitors in its class and sets a new benchmark for ride, comfort and long term reliability.

If you haven't driven it, you will not know how good it is!

Phone Holder - Bluetooth - Speakers - USB

-Charge your phone

-Play music via Bluetooth or USB stick

-Stream radio stations from your phone



Specifications*

Length	1345mm
Width	695mm
Max Load (On Flat Ground)	182kg
Max Speed up to	12km/h
Travel range up to...	40km
Suspension front & rear	YES
Tyre Sizes	300mm
Battery	45ah
Charger	6a
Motor (Taiwanese 4 pole)	
Controller PG S-Drive	120a

*subject to change without notice

Features

- * Seat Belt
- * Free Rear Basket
- * High/Low Speed Switch
- * Reclining Back Rest
- * LED Lighting & Indicators
- * Front and Rear Suspension
- * Seat Suspension