

SupaScoota® HD

(Heavy Duty)

BUILT TOUGH FOR THE HARSHTEST CONDITIONS

Finally a truly portable scooter with heaps of leg room and a large 125kg load capacity. Long distance travel in complete comfort or confined indoor use, the choice is yours. Dual motor drive system makes light work of off pavement use and steep hills so you can go anywhere with total security. A super light split frame design (heaviest part only 10kg) with seat suspension, the SupaScoota HD is designed to handle challenging conditions.

FOLDS AND LOCKS FOR TRANSPORT IN UNDER 10 SECONDS
LONG DISTANCE TRAVEL WITH OPTIONAL SPARE BATTERY PACK



Australian Designed

Unconventional and undeniably different, the SupaScoota challenges the traditional concept of portable scooter design with these exclusive features

- Acute Turning Radius
- Electronic Stability Control
- Stores in upright position
- Coil spring seat suspension
- Ultralight lithium battery option
- Twin motors for superior performance & handling

Specifications*

Length	1010mm
Width	570mm
Max Load (Level Ground)	125kg
Max Speed	6.5kph
Est. Travel Range	8 - 10km
Tyre Sizes (Solid)	200x65mm
Battery (SLA)	24v 15ah 9.2kgs
(Optional Lithium)	24v 12ah 4kgs
Motor	2 x 150w
Total Weight	29.9kg
Front Basket	Optional

*Specifications subject to change without notice

trekmobility.com.au