## PREMIUM PERFORMANCE AT AN AFFORDABLE PRICE

Extra Wide Seat - Offset Steering - 300kg Capacity

Designed from the ground up as a bespoke bariatric scooter to suit special needs clients requiring extra wide seating and a very high load rating. The *NEW* Trek Titan outperforms all competitors in it's class and sets a new benchmark for performance, ride and handling.

Optional canopy with glass windscreen provides the optimal all weather solution for golf courses and residential villages.



## Specifications\*

Length	1980mm
Width	950mm
Height	1735mm
Tyres	15in
Max Load	300kg
Max Speed programmable up to	12kph
Travel range up to	40km
Suspension front & rear	YES
Recommended Battery Size	75Ah
Charger D	C24v 8ah

\*subject to change without notice



trekmobility.com.au