



Item #	Seat	UOM
VENTURA318FS	18" Folding Seat	1/cs



HCPCS
K0807



Ventura 3 Midsize Scooter

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue
- Armrests are padded and adjustable
- Floorboard with 23" of foot space
- Ergonomic Programmable throttle control
- Easy-to-adjust delta tiller
- Height-adjustable swivel seat with fold-down backrest can be adjusted forward or backward
- Flat-free, non-marking tires and anti-tip wheels
- Rear suspension
- Headlight and large, plastic carry basket
- Easy freewheel operation

specifications			
HCPCS	K0807	Charger	4A Offboard
Weight Capacity	400 lbs.	Brakes	Electromagnetic
Top Speed	6 mph	Front Wheels	10" x 3" Flat-free
Maximum Range	16 miles	Rear Wheels	10" x 3" Flat-free
Turning Radius	53"	Base Weight	88 lbs.
Climbing Angle	7.5°	Battery Weight (pair)	28 lbs.
Ground Clearance	2.5"	Seat Weight	26 lbs.
Dimensions	48" (L) x 25" (W)	Heaviest Piece Weight	60 lbs. (Front Section)
Floor-to-Seat Height	21.5" - 23.5"	Total Weight	142 lbs.
Seat Type	Folding	Warranty on Frame	Lifetime
Seat Dimensions	18" (W) x 17.5" (D)	Warranty on Electronic Controller/Drive Train Components	24 Months
Motor	24V x 350W x 5100 rpm	Warranty on Batteries	12 Months
Batteries (pair)	12V x 33AH		