

Multifunctional electric rollator

Walker, rollator, wheelchair, mobility scooter
in one vehicle

pushing with power
assist




100KG

Mobility scooter



Pass a curb
automatically



Shopping mode with power assist

The user easily pushes the rollator loaded with heavy objects.


S
H
O
P
P
I
N
G
C
A
R
T



100KG

Automatic over curb


Hold down the power assist button , it will lift automatically and cross curb.



Product performance



Size of rollator	1000*660* (800~900)mm
Folding size	770*290*900mm
Net weight	18.6KG
Front wheel size	8 inch
Rear wheel size	8 inch
Seat size	450*290 mm
Loading weight	100KG
Frame material	Aluminum alloy
Motor power	2*250W
Speed	0~6Km/h
Cruising range	12-20Km
Battery capacity	6Ah/10AH
Charging time	2-3h
Warranty time	1 year



Power resist mode for rehabilitation exercises



Gear C1 : max resistance



Gear C2 : less resistance




Gear C3 : least resistance

Driving on the pavement, pushing it to go over the curb

Driving



Automatic over curb



Keep your Independence with...

 **SCOOTERS**
& Mobility

ABN: 13253276521

422 Stoney Crk. Rd, KINGSGROVE NSW 2208

PH: 9150 0787 FAX : 9502 3037

website: www.scootersandmobility.com.au email: dora@scootersmob.com.au