


NEW PRODUCT!



Sky Blue
Frame



Satin Grey
Frame



PRIMO 8 SEAT WALKER

RG4212BL - Sky Blue | RG4212GY - Satin Grey

FEATURES

- Great walking and stability aid
- User friendly brakes
- Height adjustable seat and handles
- Ergonomic handgrips to reduce pressure
- Comfortable padded seat
- Storage bag underneath seat.
- Padded backrest
- Simple brake adjustment
- Large diameter wheels for excellent mobility
- Flat free tyres
- Swivelling front wheels
- Folds up for transit or storage

SPECIFICATIONS

(All dimensions shown are nominal)

Overall Width: **610mm**

Internal Width: **460mm**

Depth Unfolded: **660mm**

Depth Folded: **285mm**

Overall Height: **810mm (lowest setting)**

Handle Height:
810 - 920mm (6 settings)

Seat:
360mm (W) x 350mm (D)

Seat Height:
505 - 555mm (floor to top of seat)

Wheel Diameter: **Ø200mm / 8"**

Product Weight: **8.3kg**

Max. User Weight: **136kg**



Please contact REDGUM Brand if you require further information on this product.

P: **08 9248 4180**

E: **orders@for-de.com.au**

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the For-de Group Pty Ltd

redgumbrand.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.