



FEATURES

- Great walking and stability aid
- User friendly brakes
- Height adjustable handles
- Removable basket
- Comfortable padded seat
- Padded backrest
- Large diameter wheels for excellent mobility
- Flat free tyres
- Folds up for transit or storage

SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: **700mm (W) x 700mm (D)**
770 - 920mm (H) (6 increments)
- Folded: **700mm (W) x 870mm (D)**
690mm (H)
- Seat: **510mm (W) x 360mm (D) x 540mm (H)**
- Product Weight: **11.5kg**
- Max. User Weight: **200kg**

Please contact **REDGUM Brand** if you require further information on this product



**COMFORTABLE
PADDED**
backrest and seat



Folds up for
**STORAGE
& TRANSIT**

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrand.com.au