

# **REDGUM** BRAND Daily Living and Lifestyle Products

## “CrossOver” Seat Walker & Occasional Transit Chair Walker

RG4408



### Features

- Side folding seat walker
- Converts to a transit wheelchair
- Aluminium construction
- Top quality powder coated finish
- Side folding action
- Soft operation hand brake
- Seat height at 560mm from floor
- 115mm wide padded backrest
- Removable swing away leg rests fitted
- Walking stock holder fitted as standard
- 200mm spoked wheels
- Available in three colours
- Blue / Red / /White

### Dimensions

- |                       |                 |
|-----------------------|-----------------|
| • Seat Height         | 560mm           |
| • Seat Width          | 450mm           |
| • Seat Depth          | 310mm           |
| • Unit Overall Width  | 590mm           |
| • Unit Overall Length | 980mm           |
| • Handle Height Range | 900mm to 1000mm |
| • Unit Net Weight     | 9.42kg          |

Phone: 08 9248 4180  
Fax: 08 9248 9452  
E-mail: [orders@redgumbrand.com.au](mailto:orders@redgumbrand.com.au)

**REDGUM** BRAND  
Daily Living and  
Lifestyle Products

*REDGUM Brand is a 100% Australian Owned & Operated*