

VTS SLIM DESIGN MANUAL WHEELCHAIR 863

Features:

- Compact Manual Wheelchair constructed from rust-proof 6061 aluminum alloy, offers an exceptional strength-to-weight ratio.
- Folds down in seconds with a simple pull of the seat.
- The chair features a breathable, moisture-wicking mesh seat and backrest, designed for extended wear in varying climates.
- Equipped with solid polyurethane or foam-filled tires, this wheelchair requires zero maintenance—no more flat tires on rough city sidewalks or park trails.



Technical Spec:

Material	Aluminum	Deat Width	46cm
Seat Height	49cm	Total width	64cm
Armrest	Fixed armrest	Footrest	Swing away footrest
N.W	12.3kg	G.W	14.5kg
Front wheel	6in solid tire	Rear Wheel	20in PU tire
Loading capacity	100kg	Package Size	80*30*74cm