



VOGUE WALKER RANGE

The stylish European designed Aspire Vogue Walkers are a lightweight, versatile range with a sturdy frame and easy to fold locking system, suitable for compact transportation when folded. Manufactured with the utmost attention to quality, so you can enjoy your independence with peace of mind, in style.

aspire for...

Comfort

Stability

Style



www.aspirecare.com.au

1300 133 120






DANISH DESIGN
Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Carbon - Black



Carbon - Candy Red

1
Year
Warranty

S.W.L
150kg

AS/NZ Approved

Vogue is registered in the
ARTG, AUST R ARTG 271255

MODEL

WAF705350 | Vogue Carbon Fibre - Black - Medium
WAF705360 | Vogue Carbon Fibre - Candy Red - Tall

- Luxury look and feel – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours
- Super lightweight and ergonomic design with small turning circle for excellent manoeuvrability
- AutoLocate - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate
- Locks when folded – 225mm folding width minimises storage space and makes transportation easier
- Advanced locking handbrakes – for increased user safety when walking and sitting
- Unbreakable TPE super soft wheels – offers high durability and a smooth ride for increased comfort
- Removable utility bag – simple and stylish to remove for shopping
- Supplied with cane holder attachment
- Deep Seat and Wide Backrest for improved comfort and support

	OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	WHEEL SIZE	TOTAL WEIGHT	BAG S.W.L	S.W.L
MED	695 - 920mm	680mm	600mm	225mm	240mm	460mm	550mm	8"	5kg	5kg	150kg
TALL	770 - 1040mm	680mm	605mm	225mm	240mm	460mm	620mm	8"	6kg	5kg	150kg



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Red



Champagne



S.W.L
150kg



1
Year
Warranty

Vogue is registered in the
ARTIG, AUST R ARTIG 271255

AS/NZ Approved



Deep Seat and
Wide Backrest



Cane
Holder



MODEL

- WAF705310RE | Vogue Lightweight 2 - Red
- WAF705310CH | Vogue Lightweight 2 - Champagne

- Luxury look and feel – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours
- Lightweight and ergonomic design with small turning circle for excellent manoeuvrability
- AutoLocate - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate
- Locks when folded – 225mm folding width minimises storage space and makes transportation easier
- Advanced locking handbrakes – for increased user safety when walking and sitting
- Unbreakable TPE super soft wheels – offers high durability and a smooth ride for increased comfort
- Removable utility bag – simple and stylish to remove for shopping
- Supplied with cane holder attachment
- Deep Seat and Wide Backrest for improved comfort and support

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	WHEEL SIZE	TOTAL WEIGHT	BAG S.W.L	S.W.L
785 - 920mm	655mm	620mm	225mm	240mm	465mm	540mm	8"	6.1kg	5kg	150kg



Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



White / Black



MODEL

WAF705400 | Vogue Forearm - White/Black

- Smart Posture System – equipped with forearm supports that are height, depth and angle adjustable, the user can achieve optimal centre of gravity for improved balance and weight distribution
- Rehab Focus - Easy to push with a small turning circle and low rolling resistance to accelerate standing mobility
- Locks when folded – 320mm folding width minimises storage space and makes transportation easier
- Luxury look and feel – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours
- Advanced locking handbrakes – for increased user safety when walking and sitting
- AutoLocate - Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate
- Unbreakable TPE super soft wheels – offers high durability and a smooth ride for increased comfort
- Supplied with cane holder attachment
- Removable utility bag – simple and stylish to remove for shopping

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	WHEEL SIZE	TOTAL WEIGHT	BAG S.W.L	S.W.L
1020 - 1270mm	660mm	620mm	320mm	240mm	460mm	620mm	8"	9kg	5kg	150kg




Luxury Look & Feel
Danish Designed

AVAILABLE COLOURS:



Red / Black



Champagne / Black

Super
Lightweight!
5.2kg
Meal Tray
Inclusive



Vogue is registered in the
ARTG, AUST R ARTG 271255

AS/NZ Approved



Folds to
225mm



Cane
Holder



MODEL

- WAF705100RE | Vogue Indoor - Red/Black
- WAF705100CH | Vogue Indoor - Champagne/Black

- Luxury look and feel – Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours
- Lightweight and ergonomic design with small turning circle for excellent manoeuvrability
- No lifting & spillage – meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable
- Sitting down for a meal – when the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal
- Unique folding design – the slimline design reduces the footprint when folded but allows the folded walker to stand upright on its own
- Brakes – the ergonomic ‘Soft Touch’ lockable design allows the walker to be secured in place with minimal hand pressure
- Removable storage bag – utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house
- Supplied with cane holder attachment

OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	WHEEL SIZE	TOTAL WEIGHT	BAG S.W.L	TRAY S.W.L	S.W.L
860 - 970mm	650mm	540mm	225mm	6"	5.2kg	5kg	5kg	110kg

(Single Hand Brake System and Adjustable Height IV Pole are not compatible with the Vogue Indoor Walker)



Single Hand Brake System
WAA695640



Oxygen Bottle Holder Bag
WAA695600



Night Light - LED Lamp
WAA695630



Adjustable Height IV Pole
WAA695650

NEW Aspire Vogue EXPLORER Outdoor Walker

WAF705450

COMING SOON!





ASBR21.29 220420

