



---

THE LIGHTEST WHEELCHAIR  
IN THE WORLD



# ULTRA-LITE WHEELCHAIR®

THE LIGHTEST WHEELCHAIR  
IN THE WORLD  
With top ergonomic design



## Product Specifications:

Seat Width	<b>44 cm</b>	Total Height	<b>93 cm</b>
Seat Depth	<b>41 cm</b>	Length Of Footrests	<b>39 cm</b>
Seat Height	<b>45 cm</b>	Total Length Inc. Footrests	<b>112.5 cm</b>
Backrest Height	<b>40 cm</b>	Max.User Weight	<b>120 kg</b>
Flip-up Armrest Height	<b>23.5 cm</b>	Total Weight (includes all features)	<b>11 kg</b>
Driving Unit Width	<b>68 cm</b>		